Capline

System 2010

Efficient, convenient, and versatile

The **Capline System 2010** is KNAUER's new fully automated high performance capillary electrophoresis system, designed for a wide range of different HPCE separation techniques. Combined with KNAUER's **ChromGate®** software, it is the first choice for method development tasks as well as for daily routine work.

The system contains a built-in highly sensitive UV detector to guarantee optimal performance. The wavelength is selectable at any time, even during an analysis. In addition, "stop-flow" spectra recording is supported for simplified method development.

The fully controllable liquid cooling of the capillary allows heat to dissipate more effectively and reproducibly than air-cooled systems. Thanks to an exceedingly stable voltage feed and a highly precise, pulsation-free pressure feed for hydrodynamic injections, the **Capline System 2010** offers unsurpassed reproducibility.

The 10-position inlet and outlet carousels allow for flexible handling of samples, reagents and buffers. For micro-preparative applications the outlet carousel can be used for fraction collection with KNAUER's approved fractioning algorithms.

The system is completely controlled by **ChromGate®** data system software with all of its well-known features also included. The direct control option allows for convenient and fast method development. In addition, all CE specific data are calculated automatically within the software, such as corrected peak area and mobilities.

Capline System 2010 – makes your capillary electrophoresis tasks easier than ever

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Technical Data

Voltage Range

-25 to +25 kV; maximum current: 200 µA, polarity manually reversible

Pressure Features

Highly precise, pulsation-free pressure feed with syringe pump on inlet for injection (up to 30 mbar) and rinse (1000 mbar)

Mode of Operation

Constant voltage, simultaneous freely programmable pressure

Injection Modes

Freely programmable injection and pre/post-injection tasks; hydrodynamic injection with positive pressure (up to 30 mbar); electrokinetic injection (1 to ± 25 kV)

Temperature Control of Capillary

Peltier-driven liquid cooling; temperature range: 10 °C below ambient to 50 °C; temperature stability: ± 0.1 °C

Detector

UV detector (190 - 380 nm); deuterium lamp; response time: 0.1 - 10 sec; "stop-flow" spectral scanning

Autosampler / Fraction Collector

10 position inlet carousel; 10 position outlet carousel; 1.5 ml microcentrifuge tubes; minimum sample volume: 100 µl

Safety Features

I_{min} and I_{max} programmable; safety sensors on doors to disable high voltage

Power Requirements

100 - 240V, 50 - 60 Hz, 150 W

Dimensions / Weight

420 x 360 x 440 mm (W x H x D); 24 kg

Included Software Package

ChromGate® data system with CE specific data reporting, such as time corrected peak areas and mobility. Easy operation, multi-instrument capability, FDA 21 CFR Part 11 compliance, are some of its features. Request our ChromGate® data sheet for details.

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Ordering Information

A3200	Capline System 2010, including two capillary cassettes (capillary: fused silica, ID 50 or 75 µm x 60 cm) and ChromGate® stand-alone chromatography data system
A3245	CE Start-Up Kit: vials, vial holder, syringes, buffer, ...
A1311E & A1311	KNAUER offers brand-name PCs on which ChromGate® has been completely preinstalled and tested. Our PCs meet all ChromGate® requirements and are at the current state of the art when delivered.

Technical data are subject to change without notice.

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Visit www.knauer.net for details on complete HPLC systems, HPLC columns, and osmometers.